





### Ergonomic and Versatile Design

The endoCleaner is a compact - only 47 cm wide - but ultra-technological endoscope reprocessor that meets all the safety and convenience requirements of the modern endoscopy department. The endoCleaner system allows endoscopes to be reprocessed using both ready-to-use and concentrated chemistry (HLDs). The reprocessor cleans and disinfects all types of endoscopes (gastroscopes, colonoscopes, duodenoscopes, enteroscopes, endosonographs, bronchoscopes, etc.) from various manufacturers.

- Transparent lid
- Two ergonomic control panels\*
- Integrated water preparation system
- Foot pedal for opening the lid\*
- Process documentation
- Caster wheels for easy movement

### Safety of Personnel, Patients and Equipment

Safety is one of the fundamental design principles of the reprocessor. endoCleaner has many safety options for the staff, patients, equipment and the environment.

- Automatic leak test at all stages of the reprocessing cycle
- Replaceable direct suction canisters
- Channel monitoring
- Self-disinfection

- Possibility to use single-use chemistry
- Endoscope compatibility
- Hermetically sealed lid
- UV water treatment
- Eco<sup>2</sup> System minimizes water consumption

### Chemistry Options

A unique feature of the endoCleaner Reprocessor is the dual chemistry modes: "Single Shot" mode for single-use chemistry and "Recycle" mode for multiple-use solutions.

# RFID and Data Management System

The endoCleaner can be equipped with an optional RFID system that identifies both the endoscope and the user. Network integration allows you to take advantage of the many functions of a specialized data management system.






## Fast and Advanced Treatment Cycle


The advanced cleaning system helps to significantly reduces the cleaning phase time and the concentration of the detergent required for the cycle.

#### 360° Sprinkler System



The 360° sprinkler system creates turbulence in the cleaning solution, which, through mechanical action, optimizes the cleaning of the endoscope.

#### **AntiSediment Function**




The special AntiSediment channel flushing technology prevents and eliminates blockage of the endoscope channels, guaranteeing passage for the disinfectant.

#### **Channel Monitoring Option**



The endoCleaner is optionally equipped with a signal system to monitor the flow through each channel of the endoscope, ensuring quality cleaning and disinfection during the treatment cycle.



# **Specifications**

Number of endoscopes
Cleaning method

Disinfection method

Basin capacity

Water treatment

Leak test

Alcohol flushing

Air purging

Detergents\*

Disinfectants\*

Documentation

Weight

Dimensions

Power supply

Water supply

1/2

Turbulent flushing / sprinkling / flushing

10 L - 1 endoscope / 11 L - 2 endoscopes

ilters / UV lamp

Visual / Automation

Automatic

A

Direct quotion from conjete

- Ready-to-like products: internal tank

- Concentrates: direct suction from canister Printer / internal storage / remote storage


80-100 kg

45 x 75 x 101 cm (W x H x D)

AC 220/230 V 50/60 Hz 100 W

Tap wate

### endoCleaner



\*Use only recommended chemicals. Contact your local authorized distributor for more information.







FOR COMMERCIAL AGENCIES AND MARKETING

Official Distributor: THE RHONE COMPANY FOR COMMERCIAL AGENCIES AND MARKETING

Tabasheer street, N° 12, Rawabi, Wadi Seer Amman-Jordan TEL: +962 7 7997 7494 info@rhonemed.com | www.rhonemed.com